



RNASE2 Antibody, Biotin conjugated

Product Code	CSB-PA019794LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10153
Immunogen	Recombinant Human Non-secretory ribonuclease protein (28-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-secretory ribonuclease (EC 3.1.27.5) (Eosinophil-derived neurotoxin) (RNase Upl-2) (Ribonuclease 2) (RNase 2) (Ribonuclease US), RNASE2, EDN RNS2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RNASE2