



RNASE2 Antibody, Biotin conjugated

Product Code	CSB-PA019794LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10153
Immunogen	Recombinant Human Non-secretory ribonuclease protein (28-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-secretory ribonuclease (EC 3.1.27.5) (Eosinophil-derived neurotoxin) (RNase Upl-2) (Ribonuclease 2) (RNase 2) (Ribonuclease US), RNASE2, EDN RNS2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RNASE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.