



JAM2 Antibody, HRP conjugated

Product Code	CSB-PA011936LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57087
Immunogen	Recombinant Human Junctional adhesion molecule B protein (29-238AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Junctional adhesion molecule B (JAM-B) (Junctional adhesion molecule 2) (JAM-2) (Vascular endothelial junction-associated molecule) (VE-JAM) (CD antigen CD322), JAM2, C21orf43 VEJAM
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	JAM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.