



APX2 Antibody, FITC conjugated

Product Code	CSB-PA14419C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1PER6
Immunogen	Recombinant Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic protein (4-250AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-ascorbate peroxidase 2, cytosolic (EC 1.11.1.11) (L-ascorbate peroxidase 1b) (APX1b) (AtAPx02), APX2, APX1B
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Cell Biology
Target Names	APX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.