



# Enrofloxacin Monoclonal Antibody

<b>Product Code</b>	CSB-MA00072110m
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Enrofloxacin-BSA conjugate
<b>Raised In</b>	mouse
<b>Specificity</b>	No significant cross-reactivity or interference was observed
<b>Tested Applications</b>	ELISA
<b>Form</b>	liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95% protein G purified
<b>Isotype</b>	IgG2b
<b>Clonality</b>	monoclonal
<b>Alias</b>	Baytril, enrofloxacin, CFPQ, enrofloxacin, enrofloxacinum, enrofloxacin
<b>Product Type</b>	Small Molecules Antibodies
<b>Target Names</b>	ER
<b>Accession NO.</b>	7D1F1