



Testosterone Monoclonal Antibody

Product Code	CSB-MA004411I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Testosterone-BSA conjugate
Raised In	mouse
Specificity	No significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95% ?protein G purified
Isotype	IgG2b
Clonality	monoclonal
Alias	Ts, andrusol, sustanon, teslen, testicular hormone, virosterone
Product Type	Small Molecules Antibodies
Target Names	T
Accession NO.	3E1A3