




LOXL1 Antibody

Product Code	CSB-PA013040LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08397
Immunogen	Recombinant Human Lysyl oxidase homolog 1 protein (106-574AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysyl oxidase homolog 1 (EC 1.4.3.-) (Lysyl oxidase-like protein 1) (LOL), LOXL1, LOXL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LOXL1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA013040LA01HU at dilution of 1:20

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.