



PRF1 Antibody

Product Code	CSB-PA018668LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14222
Immunogen	Recombinant Human Perforin-1 protein (22-554AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perforin-1 (P1) (Cytolysin) (Lymphocyte pore-forming protein) (PFP), PRF1, PFP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PRF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.