



KRT19 Antibody, Biotin conjugated

Product Code	CSB-PA11259D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08727
Immunogen	Recombinant Human Keratin, type I cytoskeletal 19 protein (334-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 19 (Cytokeratin-19) (CK-19) (Keratin-19) (K19), KRT19
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.