





ISG15 Antibody

Product Code	CSB-PA09744A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05161
Immunogen	Recombinant Human Ubiquitin-like protein ISG15 protein (2-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IP; Recommended dilution: WB:1:1000-1:5000, IP:1:200-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like protein ISG15 (Interferon-induced 15 kDa protein) (Interferon-induced 17 kDa protein) (IP17) (Ubiquitin cross-reactive protein) (hUCRP), ISG15, G1P2 UCRP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ISG15

Image



Western blot
 All lanes: ISG15 antibody at 2µg/ml + 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 18 kDa
 Observed band size: 18 kDa



Immunoprecipitating ISG15 in HeLa whole cell lysate

Lane 1: Rabbit control IgG instead of CSB-PA09744A0Rb in HeLa whole cell lysate. For western blotting, a HRP-conjugated Protein G antibody was used as the secondary antibody (1/5000)

Lane 2: CSB-PA09744A0Rb (8 μ g) + HeLa whole cell lysate (500 μ g)

Lane 3: HeLa whole cell lysate (20 μ g)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.