




deoD Antibody

Product Code	CSB-PA18357A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0ABP8
Immunogen	Recombinant Escherichia coli Purine nucleoside phosphorylase DeoD-type protein (2-239AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Purine nucleoside phosphorylase DeoD-type (PNP) (EC 2.4.2.1) (Inosine phosphorylase), deoD, pup
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	deoD

Image



Western blot
 All lanes: deoD antibody at 2µg/ml + DH5a whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 26 kDa
 Observed band size: 26 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.