



Cxcl14 Antibody, Biotin conjugated

Product Code	CSB-PA09319D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9WUQ5
Immunogen	Recombinant Mouse C-X-C motif chemokine 14 protein (23-99AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 14 (B-cell and monocyte-activating chemokine) (Chemokine BRAK) (Kidney-expressed chemokine CXC) (MIP-2G) (Small-inducible cytokine B14), Cxcl14, Bmac Kec Ks1 Mip2g Scyb14
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cxcl14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.