



MCAM Antibody, HRP conjugated

Product Code	CSB-PA10439B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43121
Immunogen	Recombinant Human Cell surface glycoprotein MUC18 protein (50-646AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface glycoprotein MUC18 (Cell surface glycoprotein P1H12) (Melanoma cell adhesion molecule) (Melanoma-associated antigen A32) (Melanoma-associated antigen MUC18) (S-endo 1 endothelial-associated antigen) (CD antigen CD146), MCAM, MUC18
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MCAM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.