



# Cd63 Antibody, Biotin conjugated

Product Code	CSB-PA004950ED01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41731
Immunogen	Recombinant Mouse CD63 antigen protein (103-203AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD63 antigen (CD antigen CD63), Cd63
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Cd63