



Cd63 Antibody, Biotin conjugated

Product Code	CSB-PA004950ED01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41731
Immunogen	Recombinant Mouse CD63 antigen protein (103-203AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD63 antigen (CD antigen CD63), Cd63
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Cd63
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.