



Cxcl9 Antibody, Biotin conjugated

Product Code	CSB-PA09269D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18340
Immunogen	Recombinant Mouse C-X-C motif chemokine 9 protein (22-126AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 9 (Gamma-interferon-induced monokine) (Monokine induced by interferon-gamma) (MIG) (MuMIG) (Protein m119) (Small-inducible cytokine B9), Cxcl9, Mig Scyb9
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cxcl9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.