



Cxcl9 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA09269D0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P18340 |
| Immunogen | Recombinant Mouse C-X-C motif chemokine 9 protein (22-126AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | C-X-C motif chemokine 9 (Gamma-interferon-induced monokine) (Monokine induced by interferon-gamma) (MIG) (MuMIG) (Protein m119) (Small-inducible cytokine B9), Cxcl9, Mig Scyb9 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Immunology |
| Target Names | Cxcl9 |