




IdcA Antibody, Biotin conjugated

Product Code	CSB-PA303469HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P76008
Immunogen	Recombinant Escherichia coli Murein tetrapeptide carboxypeptidase protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:2000
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Murein tetrapeptide carboxypeptidase (EC 3.4.17.13) (LD-carboxypeptidase A) (Muramoyltetrapeptide carboxypeptidase), IdcA, ycgQ
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	IdcA

Image



Western Blot

Positive WB detected in: Rosseta bacteria lysate at 40µg, 20µg, 10µg

All lanes: IdcA antibody at 1µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 34 kDa

observed band size: 34 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.