



IdcA Antibody, HRP conjugated

Product Code	CSB-PA303469HB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P76008
Immunogen	Recombinant Escherichia coli Murein tetrapeptide carboxypeptidase protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Murein tetrapeptide carboxypeptidase (EC 3.4.17.13) (LD-carboxypeptidase A) (Muramoyltetrapeptide carboxypeptidase), IdcA, ycgQ
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	IdcA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.