



Tpsab1 Antibody, Biotin conjugated

Product Code	CSB-PA024127LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02844
Immunogen	Recombinant Mouse Tryptase protein (35-267AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptase (EC 3.4.21.59) (Mast cell protease 7) (mMCP-7) (Tryptase alpha/beta-1), Tpsab1, Mcpt7
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Tpsab1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.