




Samhd1 Antibody

Product Code	CSB-PA020691LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q60710
Immunogen	Recombinant Mouse Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 protein (426-657AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxynucleoside triphosphate triphosphohydrolase SAMHD1 (dNTPase) (EC 3.1.5.-) (Interferon-gamma-inducible protein Mg11) (SAM domain and HD domain-containing protein 1), Samhd1, Mg11
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Samhd1

Image



Western blot

All lanes: Samhd1 antibody at 2µg/ml + 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 76, 75 kDa

Observed band size: 76 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.