




Ccl4 Antibody

Product Code	CSB-PA09069A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14097
Immunogen	Recombinant Mouse C-C motif chemokine 4 protein (24-92AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C motif chemokine 4 (Immune activation protein 2) (ACT-2) (ACT2) (Macrophage inflammatory protein 1-beta) (MIP-1-beta) (Protein H400) (SIS-gamma) (Small-inducible cytokine A4), Ccl4, Mip1b Scya4
Immunogen Species	Mouse
Research Area	Immunology
Target Names	Ccl4

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA09069A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.