



Alpl Antibody, Biotin conjugated

Product Code	CSB-PA17039D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09242
Immunogen	Recombinant Mouse Alkaline phosphatase, tissue-nonspecific isozyme protein (20-503AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alkaline phosphatase, tissue-nonspecific isozyme (AP-TNAP) (TNSALP) (EC 3.1.3.1) (Alkaline phosphatase 2) (Alkaline phosphatase liver/bone/kidney isozyme), Alpl, Akp-2 Akp2
Immunogen Species	Mus musculus (Mouse)
Research Area	Cell Biology
Target Names	Alpl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.