





GLA Antibody

Product Code	CSB-PA009474LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06280
Immunogen	Recombinant Human Alpha-galactosidase A protein (32-429AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-galactosidase A (EC 3.2.1.22) (Alpha-D-galactosidase A) (Alpha-D-galactoside galactohydrolase) (Melibiase) (Agalsidase), GLA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	GLA

Image



Immunohistochemistry of paraffin embedded human kidney tissue using CSB-PA009474LA01HU at dilution of 1:100



Immunohistochemistry of paraffin embedded human prostate tissue using CSB-PA009474LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.