



fliC Antibody

Product Code	CSB-PA09734A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04949
Immunogen	Recombinant Escherichia coli Flagellin protein (2-498AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Flagellin, fliC, flaF hag
Immunogen Species	Escherichia coli (strain K12)
Research Area	Immunology
Target Names	fliC