



Apoe Antibody, HRP conjugated

Product Code	CSB-PA001936LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08226
Immunogen	Recombinant Mouse Apolipoprotein E protein (19-311AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein E (Apo-E), Apoe
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Apoe
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.