



RELN Antibody, HRP conjugated

Product Code	CSB-PA019557LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78509
Immunogen	Recombinant Human Reelin protein (26-254AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reelin (EC 3.4.21.-), RELN
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RELN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.