



IpxC Antibody, HRP conjugated

Product Code	CSB-PA358910HB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A725
Immunogen	Recombinant Escherichia coli UDP-3-O-[3-hydroxymyristoyl] N-acetylglucosamine deacetylase protein (1-305AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-3-O-acyl-N-acetylglucosamine deacetylase (UDP-3-O-acyl-GlcNAc deacetylase) (EC 3.5.1.108) (Protein EnvA) (UDP-3-O-[R-3-hydroxymyristoyl]-N-acetylglucosamine deacetylase), IpxC, asmB envA
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	IpxC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.