



# VTG2 Antibody

<b>Product Code</b>	CSB-PA16369A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P02845
<b>Immunogen</b>	Recombinant Gallus gallus Vitellogenin-2 protein (26-653AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Gallus gallus
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Vitellogenin-2 (Major vitellogenin) (Vitellogenin II) [Cleaved into: Lipovitellin-1 (Lipovitellin I) (LVI); Phosvitin (PV); Lipovitellin-2 (Lipovitellin II) (LVII); YGP40], VTG2, VTGII
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Research Area</b>	Others
<b>Target Names</b>	VTG2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.