



nrdD Antibody, Biotin conjugated

Product Code	CSB-PA15989D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28903
Immunogen	Recombinant Escherichia coli Anaerobic ribonucleoside-triphosphate reductase protein (1-712AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anaerobic ribonucleoside-triphosphate reductase (EC 1.1.98.6) (Class III ribonucleoside-triphosphate reductase), nrdD
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	nrdD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.