



Anxa1 Antibody, FITC conjugated

Product Code	CSB-PA001836HC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10107
Immunogen	Recombinant Mouse Annexin A1 protein (2-346AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin I, Annexin-1, Calpactin II, Calpactin-2, Chromobindin-9, Lipocortin I, Phospholipase A2 inhibitory protein, p35, Anxa1, Anx1, Lpc-1, Lpc1, Annexin A1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Anxa1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.