



gadA Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA302526HA01ENV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P69908 |
| Immunogen | Recombinant Escherichia coli Glutamate decarboxylase alpha protein (1-466AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Glutamate decarboxylase alpha (GAD-alpha) (EC 4.1.1.15), gadA, gadS |
| Immunogen Species | Escherichia coli (strain K12) |
| Research Area | Others |
| Target Names | gadA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |