





IL6 Antibody

Product Code	CSB-PA06757A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05231
Immunogen	Recombinant Human Interleukin-6 protein (30-212AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-6 (IL-6) (B-cell stimulatory factor 2) (BSF-2) (CTL differentiation factor) (CDF) (Hybridoma growth factor) (Interferon beta-2) (IFN-beta-2), IL6, IFNB2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL6

Image



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA06757A0Rb at dilution of 1:100

**Western Blot**

Positive WB detected in: Mouse skeletal muscle tissue

All lanes: IL6 antibody at 2 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 24 kDa

Observed band size: 24 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.