




trxA Antibody

Product Code	CSB-PA18667A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AA25
Immunogen	Recombinant Escherichia coli Thioredoxin 1 protein (2-109AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin 1 (Trx-1), trxA, fipA tsnC
Immunogen Species	Escherichia coli
Research Area	Signal Transduction
Target Names	trxA

Image



Western blot

All lanes: trxA antibody at 2µg/ml + DH5α whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 12 kDa

Observed band size: 12 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.