




ligA Antibody

Product Code	CSB-PA322745HA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15042
Immunogen	Recombinant Escherichia coli DNA ligase protein (1-671AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA ligase (EC 6.5.1.2) (Polydeoxyribonucleotide synthase [NAD(+)]), ligA, dnaL lig lop pdeC
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	ligA

Image



Western blot
All lanes: DNA ligase antibody at 2µg/ml + DH5a whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 74 kDa
Observed band size: 74 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.