




Chrm1 Antibody

Product Code	CSB-PA005381HA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08482
Immunogen	Recombinant Rat Muscarinic acetylcholine receptor M1 protein (210-366AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Muscarinic acetylcholine receptor M1, Chrm1, Chrm-1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Chrm1

Image



Western Blot

Positive WB detected in: Rat brain tissue
 All lanes: Chrm1 antibody at 5.5µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 52 kDa
 Observed band size: 52 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.