




Anxa1 Antibody, Biotin conjugated

Product Code	CSB-PA001836YD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10107
Immunogen	Recombinant Mouse Annexin A1 protein (2-346AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:50-1:1000
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A1 (Annexin I) (Annexin-1) (Calpactin II) (Calpactin-2) (Chromobindin-9) (Lipocortin I) (Phospholipase A2 inhibitory protein) (p35), Anxa1, Anx1 Lpc-1 Lpc1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Anxa1

Image



Western Blot

Positive WB detected in: K562 whole cell lysate, HeLa whole cell lysate, HepG2 whole cell lysate, PC-3 whole cell lysate, Mouse lung tissue, Mouse thymus tissue

All lanes: Anxa1 antibody at 5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
predicted band size: 39 kDa
observed band size: 39 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.