



HLA-DPB1 Antibody

Product Code	CSB-PA14809A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04440
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DP beta 1 chain protein (30-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DP beta 1 chain (HLA class II histocompatibility antigen, DP(W4) beta chain) (MHC class II antigen DPB1), HLA-DPB1, HLA-DP1B
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DPB1

Image



Western blot

All lanes: HLA-DPB1 antibody at 2μg/ml

Lane 1: sw1990 whole cell lysate


Lane 2: HGC27 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA14809A0Rb at dilution of 1:50

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.