



Lox Antibody, Biotin conjugated

Product Code	CSB-PA013038LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16636
Immunogen	Recombinant Rat Protein-lysine 6-oxidas protein (163-411AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-lysine 6-oxidase (EC 1.4.3.13) (Lysyl oxidase), Lox
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Lox