



Lox Antibody, Biotin conjugated

Product Code	CSB-PA013038LD01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16636
Immunogen	Recombinant Rat Protein-lysine 6-oxidase protein (163-411AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-lysine 6-oxidase (EC 1.4.3.13) (Lysyl oxidase), Lox
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Lox
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.