



CHRNA1 Antibody, Biotin conjugated

Product Code	CSB-PA11789D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02708
Immunogen	Recombinant Human Acetylcholine receptor subunit alpha protein (21-255AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acetylcholine receptor subunit alpha, CHRNA1, ACHRA CHNRA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CHRNA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.