



CHRNA1 Antibody, Biotin conjugated

Product Code	CSB-PA11789D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02708
Immunogen	Recombinant Human Acetylcholine receptor subunit alpha protein (21-255AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acetylcholine receptor subunit alpha, CHRNA1, ACHRA CHNRA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CHRNA1