



Recombinant Human Zinc finger and BTB domain-containing protein 25 (ZBTB25)

Product Code	CSB-MP026322HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24278
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDTASHSLVL LQQLNMQREF GFLCDCTVAI GDVYFKAHRA VLAAFSNYFK MIFIHQTSEC IKIQPTDIQP DIFSYPHIM YTGKGPQKIV DHSRLEEGIR FLHADYLSHI ATEMNQVFSP ETVQSSNLYG IQISTTQKTV VKQGLEVKEA PSSNSGNRAA VQGDHPQLQL SLAIGLDDGT ADQQRACPAT QALEEHQKPP VSIKQERCDP ESVISQSHPS PSSEVTGPTF TENSVKIHL C HYCGERFDSR SNLRQHLH TH VSGSLPFGVP ASILESNDLG EVHPLNENSE ALECRRLSSF IVKENEQQPD HTNRGTTEPL QISQVSLISK DTEPVELNCN FSFSRKRKMS CTICGHKFP R KSQ LLEHMYT HKGKSYRYNR CQRFGNALAQ RFQPYCDSWS DVSLKSSRLS QEHLDLPCAL ESELTQENVD TILVE
Source	Mammalian cell
Target Names	ZBTB25
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 25 Alternative name(s): Zinc finger protein 46 Zinc finger protein KUP
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.