



Recombinant Chicken Bone sialoprotein 2 (IBSP)

Product Code	CSB-BP010945CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79780
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	FSVR SWLRRARAGD SEENAVLKSR HRYLYRYAY PPLHRYKGS SSEEGDGSE EEEEGGAPSH AGTQAAGEGL TLGDVGPGGD AASAHQDCKG GQKGTRGDSG DEDSDEEEEEE EEEEEEEEV EEQDVSVNGT STNTTAETPH GNNTVAAEEEE EDDDEEEEEE EEEEEAEAT TAAATTAQDE VTTLGDEQRS EVTTAGEQWE YEVTVGARGD EGPTESSYGD QEEPARGDSY RAYEDEYGY Y KGHGYDMYGQ DYYYNQ
Source	Baculovirus
Target Names	IBSP
Protein Names	Recommended name: Bone sialoprotein 2 Alternative name(s): Bone sialoprotein II Short name= BSP II Cell-binding sialoprotein Integrin-binding sialoprotein
Expression Region	17-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a major structural protein of the bone matrix. It constitutes approximately 12% of the noncollagenous proteins in human bone and is synthesized by skeletal-associated cell types, including hypertrophic chondrocytes, osteoblasts, osteocytes, and osteoclasts. The only extraskeletal site of its synthesis is the trophoblast. This protein binds to calcium and hydroxyapatite via its acidic amino acid clusters, and mediates cell attachment through an RGD sequence that recognizes the vitronectin receptor.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.