



Recombinant Mouse Hepatic triacylglycerol lipase (Lipc)

Product Code	CSB-MP012974MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P27656
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	<p>GGGVGTEP FGRSLGATEA SKPLKKPETR FLLFQDENDR LGCRLRPQHP ETLQECGFNS SQPLIMIIHG WSVDGLENW IWKIVSALKS RQSQPVNVL VDWISLAYQH YTIAVQNTRI VGQDVAALLL WLEESAKFSR SKVHLIGYSL GAHVSGFAGS SMDGKNKIGR ITGLDPAGPM FEGTSPNERL SPDDANFVDA IHTFTREHMG LSVGIKQPIA HYDFYPNGGS FQPGCHFLEL YKHIAEHGLN AITQTIKCAH ERSVHLFIDS LQHSDLQSIG FQCSDMGSFS QGLCLSCKKG RCNTLGYDIR KDRSGKSKRL FLITRAQSPF KVVHYQFKIQ FINQIEKPVE PTFTMSLLGT KEEIKRIPIT LGEGITSNKT YSFLITLDKD IGELILLKFK WENSAVWANV WNTVQTIMLW GIEPHHGLI LKTIWVKAGE TQQRMTFCPE NLDDLQLHPS QEKVFNCEV KSKRLTESKE QMSQETHAKK</p>
Source	Mammalian cell
Target Names	Lipc
Protein Names	Recommended name: Hepatic triacylglycerol lipase Short name= HL Short name= Hepatic lipase EC= 3.1.1.3 Alternative name(s): Lipase member C
Expression Region	23-510
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	LIPC encodes hepatic triglyceride lipase, which is expressed in liver. LIPC has the dual functions of triglyceride hydrolase and ligand/bridging factor for receptor-mediated lipoprotein uptake.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.