



ALB Antibody, Biotin conjugated

Product Code	CSB-PA001561LD01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49822
Immunogen	Recombinant Dog Serum albumin protein (25-608AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin (allergen Can f 3), ALB
Immunogen Species	Dog
Research Area	Cardiovascular
Target Names	ALB