



creN7 Antibody, Biotin conjugated

Product Code	CSB-PA502491LD01FPB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C3MPN0
Immunogen	Recombinant Sulfolobus islandicus Chromatin protein Cren7 protein (1-60AA)
Raised In	Rabbit
Species Reactivity	Sulfolobus islandicus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chromatin protein Cren7, creN7
Immunogen Species	Sulfolobus islandicus
Research Area	Others
Target Names	creN7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.