



gL Antibody, HRP conjugated

Product Code	CSB-PA515269LB01HWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F5HCH8
Immunogen	Recombinant Human cytomegalovirus Envelope glycoprotein L protein (31-278AA)
Raised In	Rabbit
Species Reactivity	Human cytomegalovirus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein L (gL), gL
Immunogen Species	Human cytomegalovirus
Research Area	Others
Target Names	gL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.