



SAA1 Antibody

Product Code	CSB-PA020656LA01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19708
Immunogen	Recombinant Dog Serum amyloid A protein (19-129AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum amyloid A protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA)], SAA1
Immunogen Species	Dog
Research Area	Others
Target Names	SAA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.