



UBE2W Antibody

Product Code	CSB-PA842638LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96B02
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 W protein (1-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 W (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme W) (N-terminal E2 ubiquitin-conjugating enzyme) (EC 2.3.2.25) (N-terminus-conjugating E2) (Ubiquitin carrier protein W) (Ubiquitin-conjugating enzyme 16) (UBC-16) (Ubiquitin-protein ligase W), UBE2W, UBC16
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBE2W
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.